

For Research Use Only

GLO1 Polyclonal antibody

Catalog Number: 18477-1-AP



Basic Information

Catalog Number:

18477-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13386

GenBank Accession Number:

BC001741

GeneID (NCBI):

2739

UNIPROT ID:

Q04760

Full Name:

glyoxalase I

Calculated MW:

21 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 18477-1-AP (GLO1 antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.