

For Research Use Only

# RASAL1 Polyclonal antibody

Catalog Number: 18504-1-AP

1 Publications



## Basic Information

Catalog Number:

18504-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13356

GenBank Accession Number:

BC014420

GeneID (NCBI):

8437

UNIPROT ID:

O95294

Full Name:

RAS protein activator like 1 (GAP1 like)

Calculated MW:

90 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Wen Wen Xu	37149929	EBioMedicine	WB, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

