

For Research Use Only

NAB1 Polyclonal antibody

Catalog Number: 18541-1-AP



Basic Information

Catalog Number:

18541-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13407

GenBank Accession Number:

BC035724

GeneID (NCBI):

4664

UNIPROT ID:

Q13506

Full Name:

NGFI-A binding protein 1 (EGR1 binding protein 1)

Calculated MW:

54 kDa

Observed MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HepG2 cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 18541-1-AP (NAB1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.