

For Research Use Only

DTX2 Polyclonal antibody

Catalog Number: 18565-1-AP **1 Publications**



Basic Information

Catalog Number: 18565-1-AP	GenBank Accession Number: BC008856	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 113878	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q86UW9	
Isotype: IgG	Full Name: deltex homolog 2 (Drosophila)	
Immunogen Catalog Number: AG4376	Calculated MW: 67 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human testis tissue, human hepatocirrhosis tissue
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xin-Yu Ding	31854042	Cancer Sci	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis using 18565-1-AP (DTX2 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 18565-1-AP (DTX2 antibody) at dilution of 1:50 (under 40x lens).