

For Research Use Only

DTX2 Polyclonal antibody

Catalog Number: 18565-1-AP

1 Publications



Basic Information

Catalog Number:

18565-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4376

GenBank Accession Number:

BC008856

GeneID (NCBI):

113878

UNIPROT ID:

Q86UW9

Full Name:

deltex homolog 2 (Drosophila)

Calculated MW:

67 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human testis tissue, human hepatocirrhosis tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xin-Yu Ding	31854042	Cancer Sci	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis using 18565-1-AP (DTX2 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 18565-1-AP (DTX2 antibody) at dilution of 1:50 (under 40x lens).