

For Research Use Only

# GADD45G-Specific Polyclonal antibody



Catalog Number: 18706-1-AP

## Basic Information

Catalog Number:

18706-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC000465

GeneID (NCBI):

10912

UNIPROT ID:

O95257

Full Name:

growth arrest and DNA-damage-inducible, gamma

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

GADD45G, also named as CR6, DDIT2 and GRP17, is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. GADD45G responds to environmental stresses by mediating activation of the p38/JNK pathway via MTK1/MEKK4 kinase. The GADD45G is highly expressed in placenta. The antibody is specific to GADD45G.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

