

For Research Use Only

GABRA6 Polyclonal antibody

Catalog Number: 18728-1-AP



Basic Information

Catalog Number:

18728-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC096241

GeneID (NCBI):

2559

UNIPROT ID:

Q16445

Full Name:

gamma-aminobutyric acid (GABA) A receptor, alpha 6

Calculated MW:

453 aa, 51 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

GABRA6 belongs to the ligand-gated ionic channel (TC 1.A.9) family. GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data