

For Research Use Only

# HOXD12 Polyclonal antibody

Catalog Number: 18735-1-AP



## Basic Information

Catalog Number:

18735-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC121103

GeneID (NCBI):

3238

UNIPROT ID:

P35452

Full Name:

homeobox D12

Calculated MW:

270 aa, 29 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue, L02 cells, mouse skeletal muscle tissue

## Background Information

HOXD12, also named as HOX4H, belongs to the Abd-B homeobox family. HOXD12 is a sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. The exact role of this gene has not been determined.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 18735-1-AP (HOXD12 antibody) at dilution of 1:200 incubated at room temperature for 1.5 hours.