

For Research Use Only

WNT8A Polyclonal antibody

Catalog Number: 18740-1-AP



Basic Information

Catalog Number: 18740-1-AP	GenBank Accession Number: BC156844	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 7478	
Source: Rabbit	UNIPROT ID: Q9H1J5	
Isotype: IgG	Full Name: wingless-type MMTV integration site family, member 8A	
	Calculated MW: 39 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse

Background Information

WNT8A, also named as WNT8D, belongs to the Wnt family of secreted signaling proteins. It is a ligand for members of the frizzled family of seven transmembrane receptors. These have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. WNT8A may play an important role in the development and differentiation of certain forebrain structures, notably the hippocampus.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data