

For Research Use Only

# ABCA8 Polyclonal antibody

Catalog Number: 18741-1-AP



## Basic Information

Catalog Number:

18741-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

bc130280

GeneID (NCBI):

10351

UNIPROT ID:

O94911

Full Name:

ATP-binding cassette, sub-family A (ABC1), member 8

Calculated MW:

1581 aa, 179 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

ABCA8 belongs to the ATP binding cassette (ABC) transporter family and it is a member of the ABC1 subfamily. ABC proteins transport various molecules across extra and intracellular membranes. The function of ABCA8 has not yet been determined.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data