

For Research Use Only

PDS5B Polyclonal antibody

Catalog Number: 18796-1-AP



Basic Information

Catalog Number:

18796-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8301

GenBank Accession Number:

BC039256

GeneID (NCBI):

23047

UNIPROT ID:

Q9NTI5

Full Name:

PDS5, regulator of cohesion maintenance, homolog B (S. cerevisiae)

Calculated MW:

1391 aa, 158 kDa, 165 kDa

Observed MW:

165-170 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HT-1080 cells, NIH/3T3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HT-1080 cells were subjected to SDS PAGE followed by western blot with 18796-1-AP (PDS5B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.