

For Research Use Only

# INTS8 Polyclonal antibody, PBS Only

Catalog Number:18802-1-PBS



## Basic Information

Catalog Number:

18802-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC136754

GeneID (NCBI):

55656

UNIPROT ID:

Q75QN2

Full Name:

integrator complex subunit 8

Calculated MW:

111 kDa

Observed MW:

95-111 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IF/ICC, IP, Indirect ELISA

Species Specificity:

human

## Background Information

Integrator complex subunit 8 (INTS8) is one of the major components of RNA polymerase II and has been demonstrated to be involved in the cleavage of small nuclear RNAs and transcriptional processes (PMID: 34880328). INTS8 was robustly increased in neurodevelopmental diseases and numerous tumors. High INTS8 expression has a tight correlation with poor prognosis in multi-cancers (PMID: 30881114). Western blot analysis detected INTS8 at an apparent molecular mass of 111 kDa.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

