

For Research Use Only

NUDT17 Polyclonal antibody

Catalog Number: 18805-1-AP



Basic Information

Catalog Number:

18805-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC046352

GeneID (NCBI):

200035

UNIPROT ID:

POC025

Full Name:

nudix (nucleoside diphosphate linked moiety X)-type motif 17

Calculated MW:

36 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

NUDT17 belongs to the Nudix hydrolase family. It is probably mediates the hydrolysis of some nucleoside diphosphate derivatives.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

