

For Research Use Only

RAB41 Polyclonal antibody

Catalog Number: 18818-1-AP

1 Publications



Basic Information

Catalog Number:

18818-1-AP

Size:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC117239

GeneID (NCBI):

347517

UNIPROT ID:

Q5JT25

Full Name:

RAB41, member RAS oncogene family

Calculated MW:

25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IP, ColP

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : PC-3 cells, human brain tissue, human kidney tissue, mouse small intestine tissue, mouse uterus tissue

Background Information

RAB41 belongs to the small GTPase superfamily and Rab family.

Notable Publications

Author	Pubmed ID	Journal	Application
Takashi Nozawa	37802980	Nat Commun	WB,IF,IP,ColP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 18818-1-AP (RAB41 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.