

For Research Use Only

RAB41 Polyclonal antibody, PBS Only

Catalog Number: 18818-1-PBS



Basic Information

Catalog Number:

18818-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC117239

GeneID (NCBI):

347517

UNIPROT ID:

Q5JT25

Full Name:

RAB41, member RAS oncogene family

Calculated MW:

25 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

RAB41 belongs to the small GTPase superfamily and Rab family.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 18818-1-AP (RAB41 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18818-1-PBS in a different storage buffer formulation.