

For Research Use Only

ADRA1D Polyclonal antibody, PBS Only

Catalog Number: 18823-1-PBS



Basic Information

Catalog Number:

18823-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000678

GeneID (NCBI):

146

UNIPROT ID:

P25100

Full Name:

adrenergic, alpha-1D-, receptor

Calculated MW:

60 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

ADRA1D, also known as ADRA1A, belongs to the G-protein-coupled receptor 1 family. ADRA1D is involved in multicellular organism aging, neuron-glia cell signaling, and regulation of blood pressure.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18823-1-AP (ADRA1D antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18823-1-PBS in a different storage buffer formulation.