

NB-3 Polyclonal antibody

Catalog Number: 18829-1-AP

Basic Information

Catalog Number:

18829-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC113118

GeneID (NCBI):

27255

UNIPROT ID:

Q9UQ52

Full Name:

contactin 6

Calculated MW:

1028 aa, 114 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

CNTN6, also named as hNB-3, belongs to the immunoglobulin superfamily and Contactin family. CNTN6 contactins mediate cell surface interactions during nervous system development. CNTN6 participates in oligodendrocytes generation by acting as a ligand of NOTCH1. Its association with NOTCH1 promotes NOTCH1 activation through the released notch intracellular domain (NICD) and subsequent translocation to the nucleus. It involved in motor coordination

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

