

For Research Use Only

GJE1 Polyclonal antibody

Catalog Number: 18851-1-AP



Basic Information

Catalog Number:

18851-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

AL033522

GeneID (NCBI):

100126572

Full Name:

gap junction protein, epsilon 1, 23kDa

Calculated MW:

24 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

GJE1, also named as Cx23, belongs to the connexin family and Beta-type (group I) subfamily. It is a connexon which is composed of a hexamer of connexins

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

