

For Research Use Only

EKLF/KLF1 Polyclonal antibody, PBS Only

Catalog Number: 18887-1-PBS



Basic Information

Catalog Number:

18887-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13467

GenBank Accession Number:

BC033580

GeneID (NCBI):

10661

UNIPROT ID:

Q13351

Full Name:

Kruppel-like factor 1 (erythroid)

Calculated MW:

362 aa, 38 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 18887-1-AP (KLF1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 18887-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using KLF1 antibody (18887-1-AP) at dilution of 1:800 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002). This data was developed using the same antibody clone with 18887-1-PBS in a different storage buffer formulation.