

For Research Use Only

HOXB9 Polyclonal antibody

Catalog Number: 18894-1-AP



Basic Information

Catalog Number:

18894-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13502

GenBank Accession Number:

BC015565

GeneID (NCBI):

3219

UNIPROT ID:

P17482

Full Name:

homeobox B9

Calculated MW:

250 aa, 28 kDa

Observed MW:

30-32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : K-562 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 18894-1-AP (HOXB9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.