

For Research Use Only

TSPAN13 Polyclonal antibody

Catalog Number: 18974-1-AP

1 Publications



Basic Information

Catalog Number:

18974-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6997

GenBank Accession Number:

BC033863

GeneID (NCBI):

27075

UNIPROT ID:

O95857

Full Name:

tetraspanin 13

Calculated MW:

204 aa, 22 kDa

Observed MW:

28-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : mouse lung tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Haofei Wang	39903772	PLoS One	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 18974-1-AP (TSPAN13 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.