

For Research Use Only

HSPA12B Polyclonal antibody, PBS Only

Catalog Number:19016-1-PBS



Basic Information

Catalog Number:

19016-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13544

GenBank Accession Number:

BC110881

GeneID (NCBI):

116835

UNIPROT ID:

Q96MM6

Full Name:

heat shock 70kD protein 12B

Calculated MW:

686 aa, 76 kDa

Observed MW:

75 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19016-1-AP (HSPA12B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 19016-1-PBS in a different storage buffer formulation.



IP result of anti-HSPA12B (IP:19016-1-AP, 4ug; Detection:19016-1-AP 1:800) with mouse brain tissue lysate 1840 ug. This data was developed using the same antibody clone with 19016-1-PBS in a different storage buffer formulation.