

For Research Use Only

FAM20B Polyclonal antibody

Catalog Number: 19134-1-AP



Basic Information

Catalog Number:

19134-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5912

GenBank Accession Number:

BC051794

GeneID (NCBI):

9917

UNIPROT ID:

O75063

Full Name:

family with sequence similarity 20, member B

Calculated MW:

46 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HepG2 cells, mouse brain tissue, mouse cerebellum tissue, rat brain tissue, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19134-1-AP (FAM20B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 19134-1-AP (FAM20B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.