

For Research Use Only

PCP4 Polyclonal antibody

Catalog Number: 19230-1-AP **1 Publications**



Basic Information

Catalog Number: 19230-1-AP	GenBank Accession Number: BC013791	Purification Method: Antigen affinity purification
Size: 250 µg/ml	GeneID (NCBI): 5121	Recommended Dilutions: IHC 1:20-1:200 IF 1:20-1:200
Source: Rabbit	UNIPROT ID: P48539	
Isotype: IgG	Full Name: Purkinje cell protein 4	
Immunogen Catalog Number: AG7247	Calculated MW: 7 kDa	

Applications

Tested Applications: IF, IHC, ELISA	Positive Controls: IHC : human cerebellum tissue, human brain tissue IF : mouse embryo,
Cited Applications: IF, IHC	
Species Specificity: human, mouse, rat	
Cited Species: mouse	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Leslie Ratié	31845734	Cereb Cortex	IHC, IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data