

For Research Use Only

# Myogenin Polyclonal antibody

Catalog Number: 19249-1-AP

6 Publications



## Basic Information

Catalog Number:

19249-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6414

GenBank Accession Number:

BC053899

GeneID (NCBI):

4656

UNIPROT ID:

P15173

Full Name:

myogenin (myogenic factor 4)

Calculated MW:

25 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human, mouse, pig

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Adrian Fischer	27901111	Sci Rep	WB
Wenxin Yin	31085189	Exp Cell Res	WB
Huiling Qiu	26624995	PLoS One	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

