

For Research Use Only

# CCDC94 Polyclonal antibody

Catalog Number: 19431-1-AP

Featured Product



## Basic Information

Catalog Number:

19431-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13764

GenBank Accession Number:

BC000561

GeneID (NCBI):

55702

UNIPROT ID:

Q9BW85

Full Name:

coiled-coil domain containing 94

Calculated MW:

323 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

IHC 1:20-1:200

## Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human, mouse, rat

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

WB : mouse brain tissue, HeLa cells

IHC : human placenta tissue,

## Background Information

CCDC94, also known as YJU2, is a splicing factor, predicted to be involved in RNA splicing and negative regulation of DNA damage response, signal transduction by p53 class mediator. Part of U2-type catalytic step 1 spliceosome. And that it is required for pre-mRNA splicing both in vivo and in vitro. Yju2 acts in concert with an unidentified heat-resistant factor(s) in an ATP-independent manner to promote the first catalytic reaction of pre-mRNA splicing after Prp2-mediated structural rearrangement of the spliceosome. (PMID: 17515604)

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




mouse brain tissue were subjected to SDS PAGE followed by western blot with 19431-1-AP (CCDC94 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human placenta using 19431-1-AP (CCDC94 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human placenta using 19431-1-AP (CCDC94 antibody) at dilution of 1:50 (under 40x lens).



WB result of CCDC94 antibody (19431-1-AP; 1:8000; incubated at room temperature for 1.5 hours) with sh-Control and sh-CCDC94 transfected HeLa cells.