

For Research Use Only

ANO9 Polyclonal antibody

Catalog Number: 19484-1-AP



Basic Information

Catalog Number:

19484-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001012302

GeneID (NCBI):

338440

UNIPROT ID:

A1A5B4

Full Name:

anoctamin 9

Calculated MW:

90 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

ANO9, also named as PIG5, TMEM16J and TP53I5, belongs to the anoctamin family. It may act as a calcium-activated chloride channel. The antibody can recognize the isoform1 and isoform3 of ANO9.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

