

For Research Use Only

NPFFR2 Polyclonal antibody, PBS Only

Catalog Number: 19505-1-PBS



Basic Information

Catalog Number:

19505-1-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_053036

GeneID (NCBI):

10886

UNIPROT ID:

Q9Y5X5

Full Name:

neuropeptide FF receptor 2

Calculated MW:

60 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Background Information

NPFFR2, also named as GPR74, NPFF2 and NPGPR, belongs to the G-protein coupled receptor 1 family. NPFFR2 is a receptor for NPAF (A-18-F-amide) and NPFF (F-8-F-amide) neuropeptides which also known as morphine-modulating peptides. NPFFr2 can also be activated by a variety of naturally occurring or synthetic FMRF-amide like ligands. NPFFR2 mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system. The antibody can recognize isoform1, isoform2 and isoform3 of NPFFR2.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human small intestine tissue slide using 19505-1-AP (NPFFR2 Antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 19505-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human small intestine tissue slide using 19505-1-AP (NPFFR2 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 19505-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 19505-1-AP (NPFFR2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 19505-1-PBS in a different storage buffer formulation.