

For Research Use Only

# IP6K3-Specific Polyclonal antibody

Catalog Number: 19518-1-AP



## Basic Information

Catalog Number:

19518-1-AP

Size:

380 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_054111

GeneID (NCBI):

117283

UNIPROT ID:

Q96PC2

Full Name:

inositol hexakisphosphate kinase 3

Calculated MW:

46 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

IP6K3, also named as IHPK3, belongs to the inositol phosphokinase (IPK) family. IP6K3 catalyze the reaction: ATP + 1D-myo-inositol hexakisphosphate = ADP + 5-diphospho-1D-myo-inositol (1,2,3,4,6)pentakisphosphate; ATP + 1D-myo-inositol 1,3,4,5,6-pentakisphosphate = ADP + diphospho-1D-myo-inositol tetrakisphosphate (isomeric configuration unknown). This antibody is specific to IP6K3

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

