

For Research Use Only

SOX17 Polyclonal antibody

Catalog Number: 19528-1-AP



Basic Information

Catalog Number:

19528-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_022454

GeneID (NCBI):

64321

UNIPROT ID:

Q9H6I2

Full Name:

SRY (sex determining region Y)-box
17

Calculated MW:

44 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

SOX17 is a transcriptional activator in the premeiotic germ cells. It binds to the sequences 5'-ACAAT-3' or 5'-AACAAAG-3'.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

