

For Research Use Only

# CDK1 Polyclonal antibody

Catalog Number: 19532-1-AP

Featured Product

220 Publications



## Basic Information

### Catalog Number:

19532-1-AP

### Concentration:

550 ug/ml

### Source:

Rabbit

### Isotype:

IgG

### GenBank Accession Number:

NM\_001786

### GeneID (NCBI):

983

### UNIPROT ID:

P06493

### Full Name:

cell division cycle 2, G1 to S and G2 to M

### Calculated MW:

34 kDa

### Observed MW:

30-34 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB 1:2000-1:10000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:100-1:400

IF/ICC 1:200-1:800

## Applications

### Tested Applications:

WB, IHC, IF/ICC, FC (Intra), IP, ELISA

### Cited Applications:

WB, IHC, IF

### Species Specificity:

human, mouse, rat, monkey

### Cited Species:

human, mouse, rat, pig, chicken, sheep

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

### Positive Controls:

**WB**: A431 cells, HeLa cells, HEK-293 cells, COS-7 cells, HepG2 cells, C2C12 cells, C6 cells, Jurkat cells, NIH/3T3 cells, Neuro-2a cells, MCF-7 cells

**IP**: HeLa cells, HEK-293 cells

**IHC**: human cervical cancer tissue, mouse liver tissue, rat liver tissue

**IF/ICC**: HeLa cells,

## Background Information

CDK1, also named as CDC2, belongs to the protein kinase superfamily, CMGC Ser/Thr protein kinase family and CDC2/CDKX subfamily. CDK1 plays a key role in the control of the eukaryotic cell cycle. It is required in higher cells for entry into S-phase and mitosis. CDK1 is a catalytic subunit of the highly conserved protein kinase complex known as M-phase promoting factor (MPF), which is essential for G1/S and G2/M phase transitions of eukaryotic cell cycle. It is a component of the kinase complex that phosphorylates the repetitive C-terminus of RNA polymerase II. Mitotic cyclins stably associate with CDK1 and function as regulatory subunits. CDK1 has 2 isoforms produced by alternative splicing with the molecular mass of 34 kDa and 27 kDa.

## Notable Publications

Author	Pubmed ID	Journal	Application
Xiaoman Chen	30207732	Mol Pharm	WB
Hui-Ting Liu	34667869	Laryngoscope Investig Otolaryngol	WB
Hongxia Li	34558311	Epigenomics	WB

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


Selected Validation Data




WB result of CDK1-Specific antibody (19532-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-CDK1-Specific transfected HeLa cells.




Various lysates were subjected to SDS PAGE followed by western blot with 19532-1-AP (CDK1-Specific antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.




Various lysates were subjected to SDS PAGE followed by western blot with 19532-1-AP (CDK1-Specific antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.




Various lysates were subjected to SDS PAGE followed by western blot with 19532-1-AP (CDK1-Specific antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.




1X10<sup>6</sup> HeLa cells were intracellularly stained with 0.4 ug Anti-Human CDK1-Specific (19532-1-AP) and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) at dilution 1:1000 (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).




Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using 19532-1-AP (CDK1-Specific antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).






Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 19532-1-AP (CDK1-Specific antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 19532-1-AP (CDK1-Specific antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CDK1-Specific antibody (19532-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 19532-1-AP (CDK1-Specific antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).

IP result of anti-CDK1-Specific (IP:19532-1-AP, 4ug; Detection:19532-1-AP 1:5000) with HeLa cells lysate 1600 ug.

IP result of anti-CDK1-Specific (IP:19532-1-AP, 4ug; Detection:19532-1-AP 1:5000) with HEK-293 cells lysate 1200 ug.