

For Research Use Only

# CD8b-Specific Polyclonal antibody

Catalog Number: 19536-1-AP

1 Publications



## Basic Information

Catalog Number:

19536-1-AP

Size:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_172213

GeneID (NCBI):

926

UNIPROT ID:

P10966

Full Name:

CD8b molecule

Calculated MW:

28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human

Cited Species:

human

## Background Information

CD8B, also named as CD8B1, identifies cytotoxic/suppressor T-cells that interact with MHC class I bearing targets. CD8 is thought to play a role in the process of T-cell mediated killing. This antibody is specific to isoform6 of CD8B.

## Notable Publications

Author	Pubmed ID	Journal	Application
Mengge Yang	36362794	J Clin Med	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data