

For Research Use Only

Phospho-PIK3R4 (Ser1079) Polyclonal antibody



Catalog Number: 19552-1-AP

Basic Information

Catalog Number:

19552-1-AP

Size:

490 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014602

GeneID (NCBI):

30849

UNIPROT ID:

Q99570

Full Name:

phosphoinositide-3-kinase,
regulatory subunit 4

Calculated MW:

153 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

PIK3R4, also named as PI3-kinase p150 subunit, Belongs to the protein kinase superfamily and Ser/Thr protein kinase family. PIK3R4 is regulatory subunit of the PI3K complex. PIK3R4 may regulate membrane trafficking late in the endocytic pathway. PIK3R4 catalyze the reaction: ATP + a protein = ADP + a phosphoprotein. The antibody is specific to Ser1079(phospho).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

