

For Research Use Only

KIAA0494 Polyclonal antibody

Catalog Number: 19569-1-AP



Basic Information

Catalog Number:

19569-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13757

GenBank Accession Number:

BC002525

GeneID (NCBI):

9813

UNIPROT ID:

O75071

Full Name:

KIAA0494

Calculated MW:

495 aa, 55 kDa

Observed MW:

45-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells, SMMC-7721 cells, HeLa cells, K-562 cells

Background Information

The function of KIAA0494 remains largely unknown.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19569-1-AP (KIAA0494 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.