

For Research Use Only

# KIAA0494 Polyclonal antibody

Catalog Number: 19569-1-AP



## Basic Information

Catalog Number:

19569-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13757

GenBank Accession Number:

BC002525

GeneID (NCBI):

9813

UNIPROT ID:

O75071

Full Name:

KIAA0494

Calculated MW:

495 aa, 55 kDa

Observed MW:

45-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : COLO 320 cells, SMMC-7721 cells, HeLa cells, K-562 cells

## Background Information

The function of KIAA0494 remains largely unknown.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19569-1-AP (KIAA0494 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.