

For Research Use Only

MCM-BP Polyclonal antibody

Catalog Number: 19573-1-AP **1 Publications**



Basic Information

Catalog Number:

19573-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13750

GenBank Accession Number:

BC004183

GeneID (NCBI):

79892

UNIPROT ID:

Q9BTE3

Full Name:

chromosome 10 open reading frame

119

Calculated MW:

640 aa, 73 kDa

Observed MW:

73 kDa, 65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuichiro Saito	35438632	Elife	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 19573-1-AP (MCM-BP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.