

For Research Use Only

# MCM-BP Polyclonal antibody

Catalog Number: 19573-1-AP

1 Publications



## Basic Information

Catalog Number:

19573-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13750

GenBank Accession Number:

BC004183

GeneID (NCBI):

79892

UNIPROT ID:

Q9BTE3

Full Name:

chromosome 10 open reading frame

119

Calculated MW:

640 aa, 73 kDa

Observed MW:

73 kDa, 65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells,

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Yuichiro Saito	35438632	Elife	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 19573-1-AP (MCM-BP antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.