

For Research Use Only

ENTPD4 Polyclonal antibody

Catalog Number: 19597-1-AP



Basic Information

Catalog Number:

19597-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5175

GenBank Accession Number:

BC034477

GeneID (NCBI):

9583

UNIPROT ID:

Q9Y227

Full Name:

ectonucleoside triphosphate
diphosphohydrolase 4

Calculated MW:

70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Jurkat cells, human placenta tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 19597-1-AP (ENTPD4 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.