

For Research Use Only

# PCDHB5 Polyclonal antibody

Catalog Number: 19609-1-AP

2 Publications



## Basic Information

Catalog Number:

19609-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7491

GenBank Accession Number:

BC001186

GeneID (NCBI):

26167

UNIPROT ID:

Q9Y5E4

Full Name:

protocadherin beta 5

Calculated MW:

86 kDa

Observed MW:

86 kDa, 70-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : HepG2 cells, mouse brain tissue, mouse liver tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yansheng Li	31367244	Theranostics	WB
Yu-Han H Hsu	37207277	iScience	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HepG2 cells lysates were subjected to SDS PAGE followed by western blot with 19609-1-AP (PCDH5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.