

# CCT6B Polyclonal antibody

Catalog Number: 19668-1-AP

## Basic Information

**Catalog Number:**

19668-1-AP

**Size:**

N/A

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_006584

**GeneID (NCBI):**

10693

**UNIPROT ID:**

Q92526

**Full Name:**

chaperonin containing TCP1, subunit 6B (zeta 2)

**Calculated MW:**

58 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human

## Background Information

CCT6B, also named as CCT-zeta-2, TCP-1-zeta-like, CCT-zeta-like, Testis-specific Tcp20 and Testis-specific protein TSA303, Belongs to the TCP-1 chaperonin family. CCT6B is a molecular chaperone which assists the folding of proteins upon ATP hydrolysis. It is known to play a role, in vitro, in the folding of actin and tubulin. The antibody has no cross reaction to CCT6A.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

