

For Research Use Only

DSG3-Specific Polyclonal antibody

Catalog Number: 19683-1-AP



Basic Information

Catalog Number:

19683-1-AP

Size:

220 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001944

GeneID (NCBI):

1830

UNIPROT ID:

P32926

Full Name:

desmoglein 3 (pemphigus vulgaris antigen)

Calculated MW:

108 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

DSG3, also named a sCDHF6 and PVA, is a component of intercellular desmosome junctions. DSG3 is involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion. This antibody is specific to DSG3.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

