

For Research Use Only

# Ephrin A4-Specific Polyclonal antibody



Catalog Number: 19685-1-AP

1 Publications

## Basic Information

Catalog Number:

19685-1-AP

Size:

170  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_182690

GeneID (NCBI):

1945

UNIPROT ID:

P52798

Full Name:

ephrin-A4

Calculated MW:

22 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

EFNA4, also named as EPLG4 and LERK4, belongs to the ephrin family. EFNA4 may play a role in the interaction between activated B lymphocytes and dendritic cells in tonsils. The antibody is specific to EFNA4.

## Notable Publications

Author	Pubmed ID	Journal	Application
Junhao Lin	34484860	Mol Ther Nucleic Acids	WB, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data