

For Research Use Only

# Jagged 2 Polyclonal antibody

Catalog Number: 19696-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

19696-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002226

GeneID (NCBI):

3714

UNIPROT ID:

Q9Y219

Full Name:

jagged 2

Calculated MW:

133 kDa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HT-29 cells, HeLa cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhikun Zheng	34421347	Int J Biol Sci	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19696-1-AP (Jagged 2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



WB result of Jagged 2 antibody (19696-1-AP; 1:400; incubated at room temperature for 1.5 hours) with sh-Control and sh-Jagged 2 transfected HeLa cells.