

For Research Use Only

# Jagged 2 Polyclonal antibody

Catalog Number: 19696-1-AP

Featured Product

2 Publications



## Basic Information

Catalog Number:

19696-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002226

GeneID (NCBI):

3714

UNIPROT ID:

Q9Y219

Full Name:

jagged 2

Calculated MW:

133 kDa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HT-29 cells, HeLa cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhikun Zheng	34421347	Int J Biol Sci	WB
Yong Feng	38485033	Gene	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19696-1-AP (Jagged 2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



WB result of Jagged 2 antibody (19696-1-AP; 1:400; incubated at room temperature for 1.5 hours) with sh-Control and sh-Jagged 2 transfected HeLa cells.