

For Research Use Only

# MYF5-Specific Polyclonal antibody

Catalog Number: 19705-1-AP

1 Publications



## Basic Information

Catalog Number:

19705-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005593

GeneID (NCBI):

4617

UNIPROT ID:

P13349

Full Name:

myogenic factor 5

Calculated MW:

28 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

pig

## Background Information

MYF5 is involved in muscle differentiation (myogenic factor). It induces fibroblasts to differentiate into myoblasts. It is probable sequence specific DNA-binding protein. The antibody is specific to MYF5.

## Notable Publications

Author	Pubmed ID	Journal	Application
Yanzhen Bi	27530319	Sci Rep	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

