

TAF4-Specific Polyclonal antibody

Catalog Number: 19725-1-AP

Basic Information

Catalog Number:

19725-1-AP

Size:

140 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003185

GeneID (NCBI):

6874

UNIPROT ID:

O00268

Full Name:

TAF4 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 135kDa

Calculated MW:

110 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

TAF4, also named as TAF2C, TAF2C1, TAF4A, TAFII130 and TAFII135, belongs to the TAF4 family. TAF4 makes part of TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. TAF4 potentiates transcriptional activation by the AF-2S of the retinoic acid, vitamin D3 and thyroid hormone. This antibody is specific to TAF4.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

