

For Research Use Only

FOXJ3 Polyclonal antibody

Catalog Number: 19751-1-AP



Basic Information

Catalog Number:

19751-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014947

GeneID (NCBI):

22887

UNIPROT ID:

Q9UPW0

Full Name:

forkhead box J3

Calculated MW:

69 kDa

Observed MW:

69 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human liver tissue, HEK-293 cells

Background Information

FOXJ3, also named as KIAA1041, contains 1 fork-head DNA-binding domain. FOXJ3 is a transcriptional activator.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 19751-1-AP (FOXJ3 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.