

For Research Use Only

BAGE2 Polyclonal antibody

Catalog Number: 19754-1-AP



Basic Information

Catalog Number: 19754-1-AP	GenBank Accession Number: NM_182482	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 85319	
Source: Rabbit	Full Name: B melanoma antigen family, member 2	
Isotype: IgG	Calculated MW: 12 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

BAGE2, also named as CT2.2, belongs to the BAGE family. BAGE2 is a candidate gene encoding tumor antigens. The antibody recognizes BAGE2 and BAGE3.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data