

OR10J5 Polyclonal antibody

Catalog Number: 19761-1-AP

Basic Information

Catalog Number:

19761-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001004469

GeneID (NCBI):

127385

UNIPROT ID:

Q8NHC4

Full Name:

olfactory receptor, family 10,
subfamily J, member 5

Calculated MW:

34 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

OR10J5, also named as OR1-28, belongs to the G-protein coupled receptor 1 family. OR10J5 is an odorant receptor.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

