

For Research Use Only

POLQ Polyclonal antibody

Catalog Number: 19770-1-AP



Basic Information

Catalog Number:

19770-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_199420

GeneID (NCBI):

10721

UNIPROT ID:

O75417

Full Name:

polymerase (DNA directed), theta

Calculated MW:

124 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

POLQ, Belongs to the DNA polymerase type-A family. POLQ could be involved in the repair of interstrand cross-links. POLQ catalyze the reaction: Deoxynucleoside triphosphate + DNA(n) = diphosphate + DNA(n+1).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

